

BULLETIN DETAILS

Product Group: Quantum Powerchairs

Effective Date:

22 September 2025

Model(s) Affected: Edge 3 Stretto

Description of Changes/ Diviations:

This bulletin informs providers of a design update to the battery door panel on the Edge 3 Stretto.

Details of Change

- Old Style:
 - Battery door panel with an opening mechanism located on the front (see image).
 - Panel colour matched to shroud
- New Style:
 - Redesigned battery door panel with two side screws securing the door.
 - Door panel will now be supplied in black only and will no longer be colour-matched to the shroud
 - The updated design ensures the door remains more secure and stable during use.

OLD BATTERY DOOR



OLD OPENING MECHANISM CIRCLED

NEW BATTERY DOOR



NEW OPENING MECHANISM CIRCLED

Documentation Affected: Owner's Manual Order Form(s) Spec Sheets Other N/A

Link to new version of affected document:

For more detailed information, please contact your local designated Area Sales Manager.

! NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

Pride Mobility Products Ltd.
32 Wedgwood Road, Bicester,
Oxfordshire, OX26 4UL

E: sales@quantumrehab.co.uk
T: +44 (0)1869 324600