



NEW Smart Traction Control (STC)



Q-Logic 3e Controller



Superior, Automotive-Grade SRS (Smooth Ride Suspension)



LED front and rear fender lights

4FRONT²

The 4Front[®] 2 is the next generation of the 4Front Power Chair! New Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), provide stability for smooth obstacle transitions and unmatched driving performance. Additional new features include 4-pole motors and redesigned caster arms for improved durability and aesthetics and reduced weight.

01869 324600 | quantumrehab.co.uk



Specifications

Power Chair

Weight Capacity ⁴	21 st 6 lbs (136 kg)
Maximum Speed ¹	Up to 6.25 mph / Up to 3.5 mph when Elevated
Ground Clearance	80 mm (3.15") at Battery Tray
Turning Radius ²	657.86 mm (25.9")
Overall Size ²	Length: 1013.46 mm (39.9") without Front Riggings Width: 609.6 mm (24")
Component Weights ⁷	Base: 130 lbs (58.9 kg) Batteries: 51.5 lbs (23.36 kg) / 44.8 lbs (20.32 kg)
Available Eletronics ⁸	120A Q-Logic 3e (Standard) 120A Q-Logic 3 (Optional)
Moter Packages	4-Pole, Front-Wheel Drive

Wheels

Drive Wheels	Standard: 355.6 mm (14") Flat Free Optional: 355.6 mm (14") Knobby/Pneumatic
Caster Wheels	Rear: 228.6 mm (9") Flat Free
Anti-Tip Wheels	Front: 101.6 mm (4") Flat Free
Suspension	Superior, Automotive-Grade SRS (Smooth Ride Suspension)

Batteries

Batteries ³	70Ah Batteries
Battery Charger	8-Amp, Off-Board

Seating Sizes

TRU-Balance® 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Lift & Tilt	Width: 12-22" Depth: 12-22"

Seat to Floor Heights

TRU-Balance® 3 Static Seat ⁵	15.2-17.25"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Power Lift & Tilt ⁵	17.75-19.75"

Available Options

10" of Power Adjustable Safe Seat Elevation at Speeds of up to 3.5 mph

NEW Backup Camera

Trilogy Vent Tray ⁸

Compact Vent Tray ⁸

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings.

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Battery weight may vary +/- 3 lbs. based on manufacturer.

7) All Quantum Rehab® units can be ordered with specialty controls.

8) Call for availability.

9) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front 2's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

Colours

Please contact your representative for the full range of options and colours.



Blue
Jean



Red
Corvette



Lemon
Crush



Green
Machine



Raspberry
Beret



Evolution
Orange



Sugar
Plum



REV 01 - 25/07/22



quantumrehab.co.uk
01869 324600