



# QUANTUM 1450

The Quantum® 1450 front-wheel drive bariatric power base features innovative design that provides exceptional outdoor performance while delivering outstanding tight space maneuverability. The 1450 accepts a complete range of rehab seating and electronics options.

QUANTUM®

## Specifications

<b>WEIGHT CAPACITY</b> <sup>4</sup>	42 stone 9 lbs. (270 kg)
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 4 mph
<b>GROUND CLEARANCE</b> <sup>2</sup>	38 mm (1.5") at anti-tip wheels
<b>TURNING RADIUS</b> <sup>2</sup>	673 mm (26.5") with front riggings
<b>OVERALL SIZE</b> <sup>2, 6</sup>	<b>Length:</b> 971 mm (38.25") without front riggings 1149 mm (45.25") with front riggings <b>Width:</b> 743 mm (29.25")
<b>DRIVE WHEELS</b>	356 mm (14") solid
<b>CASTER WHEELS</b>	<b>Front:</b> 229 mm (9") solid <b>Rear:</b> 229 mm (9") solid
<b>SUSPENSION</b>	Swing Arm Suspension
<b>COMPONENT WEIGHTS</b> <sup>5, 7</sup>	<b>Base:</b> 195 lbs (88.5 kg) <b>Batteries:</b> 52 lbs (23.6 kg) each
<b>DRIVETRAIN</b>	Two motor, front wheel drive
<b>BRAKES</b>	"Intelligent Braking" electronic regenerative, disc park brake
<b>BATTERIES</b> <sup>3</sup>	Two 12-volt, deep-cycle, 70 Ah batteries
<b>BATTERY CHARGER</b>	8-amp, Off-board
<b>MAXIMUM OBSTACLE CLIMBING ABILITY</b>	60 mm (2.36")
<b>SEATING SIZES</b>	<b>TRU-Balance® 2 HD Tilt</b>
<b>Synergy®/Static Seat</b>	<b>Width:</b> 20-32" <b>Depth:</b> 20-28"
<b>Power Tilt</b>	<b>Width:</b> 20-32" <b>Depth:</b> 20-28"
<b>Lift &amp; Tilt</b>	N/A
<b>SEAT-TO-FLOOR HEIGHTS</b>	<b>TRU-Balance® 2 HD Tilt</b>
<b>Synergy®/Static Seat</b> <sup>8</sup>	16.5-18.5"
<b>Power Tilt</b> <sup>8</sup>	19.5-20.5"
<b>Lift &amp; Tilt</b> <sup>8</sup>	N/A

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
 Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability


## Features

- 42 stone 9 lbs. weight capacity
- Maximum speed up to 4 mph
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems

## Options

- TRU-Balance® 2 Power Tilt
- Manual and Power Elevating Leg Rests
- Comfort seating
- VR2 Controls
- NE+ Controls
- Q-Logic 3 Controls

## Colours



\* Please contact your representative for the full range of options and colours.



**QUANTUM®**

01869 324600

sales@quantumrehab.co.uk

quantumrehab.co.uk

