

# Technical Bulletin

Bulletin Ref.:

UK-1-2-26

Date Issued:

10 February 2026

## BULLETIN DETAILS

**Product Group:** Quantum Powerchairs

**Effective Date:**

10 February 2026

**Model(s) Affected:** Quantum Outback

### Description of Changes/ Diviations:

This bulletin provides information regarding a charger update for the Quantum Outback powerchair. The existing 8 amp charger is being replaced with a 10 amp charger to support improved charging performance. This change will be effective from 10 February 2026. All Quantum Outback powerchairs manufactured from this date onwards will be supplied with the 10 amp charger.

### Service and Support

- Replacement chargers should be ordered as the 10 amp version from the effective date
- No changes to charging procedure are required

### Dealer / Technician Action

- Update parts ordering and service records to reflect the new charger specification
- Inform service and workshop teams of the change



**OLD 8 AMP CHARGER**



**NEW 10 AMP CHARGER**

**Documentation Affected:**  Owner's Manual  Order Form(s)  Spec Sheets  Other  N/A

**Link to new version of affected document:**

*For more detailed information, please contact your local designated Area Sales Manager.*

NB: This Bulletin is for informational purposes only and DOES NOT replace the Specifications.

**Pride Mobility Products Ltd.**  
32 Wedgwood Road, Bicester,  
Oxfordshire, OX26 4UL

E: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)  
T: +44 (0)1869 324600